

CERTIFICATE OF ANALYSIS

Proteinase K from Tritirachium album lyophil.	Cat.No. : 33752 Contr.No.: 150294
---	--

Parameter	Method	Specification	Result
Molecular weight		ca. 28 000	
Appearance		white lyophilisate	corresponds
Activity	DMC U/mg determined with TNBS	min. 8	11.0
Extraneous activities	DNase RNase	not detectable not detectable	corresponds corresponds
Minimum shelf life			03/ 2017
Storage (°C)			+2 to +8

Unit definition

One DMC-Unit is that amount of enzymatic activity which catalyzes the cleavage of 1 µequivalent peptide bond from dimethylcasein per minute at 25°C, pH 7.0, expressed in terms of the appearance of new terminal amino groups.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 25.03.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.