

CERTIFICATE OF ANALYSIS

Proteinase K from Tritirachium album
lyophil.

Cat.No. : 33752
Contr.No.: 230006

Parameter	Method	Specification	Result
Molecular weight		ca. 28 000	
Appearance		white lyophilisate	corresponds
Activity (mAnsonU/mg)		min. 30	> 30
Extraneous activities	DNase	not detectable	corresponds
	RNase	not detectable	corresponds
Minimum shelf life			30.04.2025
Storage (°C)			+2 to +8

Unit definition

One unit is defined as the amount of enzyme that liberates Folin-positive amino acids and peptides, corresponding to 1 pmole tyrosine under assay conditions in 1 minute using hemoglobin as substrate.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 11.03.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.