

CERTIFICATE OF ANALYSIS

Proteasome Inhibitor MG-132	Cat. No. : 33766
	Contr. No. : 140086

Parameter	Method	Specification	Result
Amino acid-sequence		Z-Leu-Leu-Leu-CHO	
Formula		C ₂₆ H ₄₁ N ₃ O ₅	
Molecular weight		475.62	
Appearance		white lyophilisate	corresponds
Assay (%)	HPLC	min. 90.0	corresponds
Purity (%)	from N elemental analysis		99.5
Solubility		soluble in DMSO	corresponds
Advise		do not use an aqueous buffer substitute for DMSO because this compound may racemize soon	
Minimum shelf life			12.01.2018
Storage / Stability	powder	at -15 to -25 °C in the dark for more than one year	
	solution	at -15 to -25 °C for more than one month if it is protected from moisture and kept in the dark	

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

**SERVA Electrophoresis GmbH
Quality Control**

Printing date: 17.10.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.