

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Detergent 7x neutral, phosphatfrei	Kat.Nr./Cat.No. : 34205
neutral, phosphate-free	Lot/Contr.No. : P110467

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weiß-gelbliche, viskose Flüssigkeit white-yellowish, viscous liquid	entspricht corresponds
Anwendung use		Reinigung cleaning procedures	
pH	1:10 verd./ dil. in H ₂ O	7.0 – 8.0	7.3
Leitfähigkeit (mS/ cm) conductivity	1:10 verd./ dil. in H ₂ O	max. 5.0	3.2
Mindesthaltbarkeit minimum shelf life			11/ 2016
Lagerung (°C) storage			+15 bis/ to +30

Farbunterschiede von lot zu lot sind möglich.
The colour may change from batch to batch.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.11.2011

Dipl.-Ing. (FH) Bernhard Göckel

Dipl.-Ing. (FH) Christian Nowak

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.