

CERTIFICATE OF ANALYSIS

| | |
|--|---------------------------|
| Ribonuclease A from bovine pancreas | Cat.No. : 34388 |
| E.C. 3.1.27.5 lyophil. | Contr.No. : 131298 |

| Parameter | Method | Specification | Result |
|---------------------------|-----------------------------|----------------------|---------------|
| Molecular weight | | ca. 13 700 | |
| Appearance | | white lyophilisate | corresponds |
| Solubility | 1 mg/ ml H ₂ O | clear and colorless | corresponds |
| Activity (U/ mg) | Ribonuclease A by Kunitz | min. 80 | 91 |
| DNase activity | | none detected | corresponds |
| Minimum shelf life | | | 07/ 2015 |
| Storage (°C) | | | +2 to +8 |

Unit definition

1 Kunitz-Unit is that amount of enzymatic activity which is capable of causing within 1 minute a decrease in absorbance at 300 nm equivalent to the maximum possible change in a 0.05 % solution of yeast RNA at 25 °C, pH 5.0.

| |
|--|
| We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods. |
|--|

SERVA Electrophoresis GmbH
Quality Control**Printing date: 25.11.2013**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.