

**CERTIFICATE OF ANALYSIS**

<b>Ribonuclease A from bovine pancreas</b> lyophil.	<b>Cat.No. : 34390</b> <b>Contr.No.: 190534</b>
--	--

Parameter	Method	Specification	Result
<b>Molecular weight</b>		ca. 13 700	
<b>Appearance</b>		white lyophilisate	corresponds
<b>Solubility</b>	1 mg/ml H <sub>2</sub> O	clear and colorless	corresponds
<b>Activity (U/mg)</b>	Ribonuclease A by Kunitz	min. 70	152
<b>Minimum shelf life</b>			28.02.2023
<b>Storage (°C)</b>			+2 to +8

**Unit definition**

1 Kunitz-Unit is that amount of enzymatic activity which is capable of causing within 1 minute a decrease in absorbance at 300 nm equivalent to the maximum possible change in a 0.05 % solution of yeast RNA at 25 °C, pH 5.0.

**We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.**

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date: 24.05.2019**

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.