

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS**SERVA Violet 17 Staining Kit****Kat.Nr./Cat.No. : 35074****Lot/Contr.No. : P110532**

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	zum Färben von IEF-Gelen for staining of IEF gels	
<b>Inhalt</b> content	Violet 17 tablets	Cat.no. 35075
	20% TCA Solution for Fixation	Cat.no. 36913
	Stain Solubilizer Solution (20% phosphoric acid)	Cat.no. 42562
	Destaining Solution (3% phosphoric acid)	Cat.no. 42563
<b>Mindesthaltbarkeit</b> minimum shelf life		02/ 2014
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control****Druckdatum/ printing date: 07.02.2013**

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.