

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVA Violet 17 Staining Kit

Kat.Nr./Cat.No. : 35074

Lot/Contr.No. : P130368

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	zum Färben von IEF-Gelen for staining of IEF gels	entspricht corresponds
Inhalt content	Violet 17 tablets 20% TCA Solution for Fixation Stain Solubilizer Solution (20% phosphoric acid) Destaining Solution (3% phosphoric acid)	
Mindesthaltbarkeit minimum shelf life		02/ 2015
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.09.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.