

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVA Violet 17 Staining Kit

Kat.Nr./Cat.No. : 35074

Lot/Contr.No. : P150135

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	zum Färben von IEF-Gelen for staining of IEF gels	entspricht corresponds
Inhalt content	1 g SERVA Violet 17 500 ml 20% TCA Solution for Fixation 500 ml Stain Solubilizer Solution (20% phosphoric acid) 2 x 500 ml Destaining Solution (3% phosphoric acid)	
Mindesthaltbarkeit minimum shelf life		10/ 2016
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 23.04.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.