

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

SERVA Violet 17 Staining Kit	Kat.Nr./Cat.No. :	35074
	Lot/Contr.No. :	P170097

Parameter parameter	Spezifikation specification	Ergebnis result
Anwendung application	zum Färben von IEF-Gelen for staining of IEF gels	entspricht corresponds
Inhalt content	1 g SERVA Violet 17 500 ml 20% TCA Solution for Fixation 500 ml Stain Solubilizer Solution (20% phosphoric acid) 2 x 500 ml Destaining Solution (3% phosphoric acid)	
Mindesthaltbarkeit minimum shelf life		06.10.2018
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 06.03.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.