

**ANALYSENZERTIFIKAT**  
CERTIFICATE OF ANALYSIS

---

<b>SERVA Violet 17 Staining Kit</b>	<b>Kat.Nr./Cat.No. :</b>	<b>35074</b>
	<b>Lot/Contr.No. :</b>	<b>P170297</b>

---

<b>Parameter</b> parameter	<b>Spezifikation</b> specification	<b>Ergebnis</b> result
<b>Anwendung</b> application	zum Färben von IEF-Gelen for staining of IEF gels	entspricht corresponds
<b>Inhalt</b> content	1 g SERVA Violet 17 500 ml 20% TCA Solution for Fixation 500 ml Stain Solubilizer Solution (20% phosphoric acid) 2 x 500 ml Destaining Solution (3% phosphoric acid)	
<b>Mindesthaltbarkeit</b> minimum shelf life		01.12.2018
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

---

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 06.09.2017**

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.