

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

| | |
|---|--|
| SERVA Silver Stain Kit for SDS PAGE | Kat.Nr./Cat.No. : 35076 Lot/Contr.No. : P160092 |
|---|--|

| Parameter parameter | Spezifikation specification | Ergebnis result |
|---|--|---------------------------|
| Anwendung application | SDS PAGE | entspricht corresponds |
| Inhalt content | 500 ml Fixing solution 25 x 30mg Natriumthiosulfate 5 H ₂ O 250 ml Silver nitrate solution 1.5 ml Formaldehyde 500 ml Sodium carbonate solution 30 g Glycine | |
| Mindesthaltbarkeit minimum shelf life | | 22.06.2017 |
| Lagerung (°C) storage | | +15 bis/ to +30 |

**Bei Lagerung unter +10 °C kann die Natriumcarbonat Lösung auskristallisieren.
Die Kristalle lösen sich beim Erwärmen des Produktes wieder.**

Storage below +10 °C may cause crystallization of Sodium carbonate solution.
Crystals can be solved by warming up the product.

| |
|--|
| <p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p> |
|--|

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 08.03.2016

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.