

ANALYSENZERTIFIKAT

CERTIFICATE OF ANALYSIS

SERVA Silver Stain Kit
for Native PAGE

Kat.Nr./Cat.No. : 35077
Lot/Contr.No. : P220421

Parameter
parameter

Spezifikation
specification

Ergebnis
result

Anwendung
application

native PAGE

Inhalt
content

500 ml Trichloroacetic acid
25 x 30mg Sodiumthiosulfate· 5 H₂O
250 ml Silver nitrate solution
1.5 ml Formaldehyde
500 ml Sodium carbonate solution
30 g Glycine

Mindesthaltbarkeit
minimum shelf life

30.09.2023

Lagerung (°C)
storage

+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 02.12.2022

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.