

A N A L Y S E N Z E R T I F I K A T
CERTIFICATE OF ANALYSIS

SERVA ProteinStain Fluo-Y, Lösung 100 x	Kat.Nr./Cat.No. : 35092
SERVA ProteinStain Fluo-Y, solution, 100 x	Lot/Contr.No. : 130931

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		gelbe, viskose Lösung yellow, viscous solution	entspricht corresponds
Anwendung application	Proteindetektion in SDS-PAGE protein detection in SDS-PAGE	entspricht corresponds	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			14.08.2017
Lagerung (°C) storage			+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 20.06.2017

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.