

A N A L Y S E N Z E R T I F I K A T
 CERTIFICATE OF ANALYSIS

Silikon Anti-Foam Emulsion, 30 %, USP	Kat.Nr./Cat.No. :	35119
Silicone anti-foam emulsion, 30 %, USP	Lot/Contr.No. :	160959

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		weiße Emulsion white emulsion	entspricht corresponds
Gehalt (%) assay	Polydimethylsiloxane	28.0 – 32.0	29.5
pH		2.0 – 3.0	2.8
Antischaumaktivität (s) defoaming activity	USP	max. 15	3
Schwermetalle (ppm) heavy metals	USP	max. 5	entspricht corresponds
Mindesthaltbarkeit minimum shelf life			10.05.2019
Lagerung (°C) storage			+15 bis / to +30

Dieses Produkt enthält keine Bestandteile tierischen Ursprungs.
 This product does not contain any components of animal origin.

Die physikalisch-chemischen Eigenschaften entsprechen den Anforderungen der USP.
 The physical/ chemical data correspond to the requirements of USP.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
 We do not guarantee that the product can be used for a special application.
 This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 25.11.2016

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.