

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Siliconlösung SERVA in Isopropanol	Kat.Nr./Cat.No. : 35130
Silicone solution SERVA in isopropanol	Lot/Contr.No. : P180191

Parameter parameter	Spezifikation specification	Ergebnis result
Aussehen appearance	klare, farblose Flüssigkeit clear and colorless liquid	entspricht corresponds
Anwendung use	Zur Siliconisierung von Glas und Metall geeignet. Tested for siliconizing of glass and metal.	
Mindesthaltbarkeit minimum shelf life		05.06.2023
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 07.06.2018

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.