

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Siliconöl DC 200; 10 cst pract.	Kat.Nr./Cat.No. : 35132
Silicone DC 200 fluid; 10 cst, pract.	Lot/Contr.No. : 170591

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, farblose Flüssigkeit clear, colorless liquid	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Dichte (g/ml) specific gravity	1013 hPa/ 25 °C	0.934 - 0.940	entspricht corresponds
Brechungsindex refractive index	25 °C	1.3980 - 1.4000	1.3990
Viskosität (cSt) viscosity		9.0 - 11.0	9.6
Flüchtige Bestandteile (%) volatile content		max. 1.0	0.16
Mindesthaltbarkeit minimum shelf life			08.12.2020
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 18.11.2019

Christian Monsler

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.