

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Siliconöl DC 200; 100 cst pract.	Kat.Nr./Cat.No. :	35135
Silicone DC 200 fluid; 100 cst pract.	Lot/Contr.No. :	120172

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, viskose Flüssigkeit clear, viscous liquid	entspricht corresponds
Identität identity	IR	entspricht corresponds	entspricht corresponds
Dichte (g/ml) specific gravity		0.962 – 0.968	entspricht corresponds
Brechungsindex refractive index	25 °C	1.4020 – 1.4040	1.4027
Viskosität (cSt; 25 °C) viscosity		95.0 – 105.0	99.5
Flüchtige Bestandteile (%) volatile content		max. 0.5	≤ 0.5
Mindesthaltbarkeit minimum shelf life			12/ 2014
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.02.2012

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.