

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Siliconöl DC 710; 500 cst pract.	Kat.Nr./Cat.No. :	35149
Silicone DC 710 fluid; 500 cst pract.	Lot/Contr.No. :	230394

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Aussehen appearance		klare, viskose Flüssigkeit clear, colorless liquid	entspricht corresponds
Dichte (g/ml) specific gravity	25 °C	1.10 – 1.11	entspricht corresponds
Brechungsindex refractive index	25 °C	1.5300 – 1.5380	1.5359
Viskosität (cst) viscosity	25 °C	450.0 – 550.0	535
Mindesthaltbarkeit minimum shelf life			26.03.2027
Lagerung (°C) storage			+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 12.04.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.