

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Synperonic® F108 pract. **Kat.Nr./Cat.No. : 35726**
Lot/Contr.No. : 140475

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Molekulargewicht molecular weight		ca. 14 000	
Aussehen appearance		weiße Schuppen white flakes	entspricht corresponds
Löslichkeit solubility	2.5 % in H ₂ O	klare Lösung clear solution	entspricht corresponds
Farbe (Hazen) colour	25 % in H ₂ O	max. 50	12
Schmelzpunkt (°C) melting point		55 – 60	59
Hydroxylzahl (mg KOH/g) hydroxyl number		6.5 – 8.9	8.4
pH	2.5 % in H ₂ O	5.0 – 7.5	7.2
Wasser (%) water	KF	max. 0.75	0.15
Sulfatasche (%) sulfated ash		max. 0.2	0.1
Mindesthaltbarkeit minimum shelf life			12/ 2015
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.
We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 15.04.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
This report has been computer generated and does not contain a signature.