

ANALYSENZERTIFIKAT
CERTIFICATE OF ANALYSIS

Synperonic® F108 pract.

Kat.Nr./Cat.No. : 35726

Lot/Contr.No. : 240066

Parameter parameter	Methode method	Spezifikation specification	Ergebnis result
Molekulargewicht molecular weight		ca. 14 000	
Aussehen appearance		weiße Schuppen white flakes	entspricht corresponds
Löslichkeit solubility	2.5 % in H ₂ O	klare Lösung clear solution	entspricht corresponds
Farbe (Hazen) colour	25 % in H ₂ O	max. 50	16
Schmelzpunkt (°C) melting point		55 – 60	58.8
Hydroxylzahl (mg KOH/g) hydroxyl number		6.5 – 8.9	8.2
pH	2.5 % in H ₂ O	5.0 – 7.5	7.0
Wasser (%) water	KF	max. 0.75	0.28
Sulfatasche (%) sulfated ash		max. 0.2	0.1
Mindesthaltbarkeit minimum shelf life			12.12.2024
Lagerung (°C) storage			+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 31.01.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.