

# CERTIFICATE OF ANALYSIS

<b>Thrombin from bovine plasma</b> min. 1000 units/mg protein lyophil. E.C. 3.4.21.5	<b>Cat. No. :</b> 36402 <b>Lot. No. :</b> 200376
---	---

Parameter	Method	Specification	Result
<b>Molecular weight</b>		ca. 37 000	
<b>Appearance</b>		white lyophilisate	corresponds
<b>Specific Activity (NIH U/mg protein)</b>	Thrombin	min. 1000	2659
<b>Activity (U/mg solid)</b>	Thrombin		264.42
<b>Minimum shelf life</b>			07.04.2023
<b>Storage (°C)</b>			+2 to +8

### Unit definition:

One NIH unit clots a standard fibrinogen solution in 15 sec. at 37 °C.

This material is derived from animals born and sourced from the United States. The animals are from a natural beef program in accordance to FDA Regulation 21 CFR 589.2000, which prohibits the use of animal proteins in ruminant feed.

The animals, both pre- and post-mortem, were apparently free from infectious and contagious diseases. Therefore to the best of our knowledge and belief, we can state that BSE/TSE should not be a concern.

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b>  <b>Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.</b>          We do not guarantee that the product can be used for a special application.          This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date:** 08.04.2020

Christian Monsler

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
 This report has been computer generated and does not contain a signature.