

CERTIFICATE OF ANALYSIS

Thrombin from bovine plasma	Cat. No. : 36402
min. 1000 units/mg protein lyophil. E.C. 3.4.21.5	Lot. No. : 240288

Parameter	Method	Specification	Result
Molecular weight		ca. 37 000	
Appearance		white lyophilisate	corresponds
Specific Activity (NIH U/mg protein)	Thrombin	min. 1000	2019
Activity (U/mg solid)	Thrombin		601
Origin			New Zealand
Minimum shelf life			21.03.27
Storage (°C)			+2 to +8

Unit definition:

One NIH unit clots a standard fibrinogen solution in 15 sec. at 37 °C.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 22.03.2024

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.