

CERTIFICATE OF ANALYSIS

PNGase F, recombinant	Cat.No.: : 36404
	Contr.No.: 210086

Parameter	Method	Specification	Result
Molecular weight	PNGase F	ca. 36 000	
Purity (%)	SDS-PAGE	≥ 95	≥ 95
Activity (U/ml)		1 000 000	1 000 000
Form	The enzyme is supplied as a liquid in 1 x PBS (Phosphate Buffered Saline; 137 mM NaCl, 10 mM Phosphate, 2.7 mM KCl, pH 7.4)		
Origin		E. coli	corresponds
Minimum shelf life			27.01.2023
Storage (°C)			-15 bis/ to -25

Unit definition:

Achieves complete deglycosylation of 10 µg of RNase B incubated in 1 x PBS with 1 µl of PNGase F for 5 – 10 minutes at 37 °C or room temperature. Separation of reaction products are visualized by SDS-PAGE.

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 02.02.2021

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.