

CERTIFICATE OF ANALYSIS

Transferrin human (Apo) lyophil.	Cat.No. : 36760 Contr.No.: 090381
--	--

Parameter	Method	Specification	Result
Molecular weight		ca. 77 000	
Appearance		off-white to slightly pink lyophil.	corresponds
Protein content (%)	calc.on dried basis	min. 95,0	97.0
pH	1% solution	6.5 – 8.0	6.8
Total bound Iron (µg/g)	Atomic Absorption	max. 30	30
Total Iron Binding Capacity (mg/g protein)		min. 1.0	2.0
IgG – free	at 2 %	negative for anti-IgG	negative
Endotoxin (EU/mg)	Limulus Amebocyte Lysate	max. 1.0	< 1.0
Moisture (%)	Loss on drying	max. 3.0	0.4
Minimum shelf life			10/ 2014
Storage (°C)			-15 to -25

This product was heat treated at 60°C for 10 hours for virus inactivation prior to lyophilization. Each unit of plasma used in the manufacture of this product has been tested and found to be non-reactive for Hepatitis B Surface Antigen (HbsAg), antibody to HIV I, II and anti-HCV, and Syphilis by FDA approved tests.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 05.08.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.