

CERTIFICATE OF ANALYSIS

Transferrin human (Apo) lyophil.	Cat.No. : 36760 Contr.No.: 180648
--	--

Parameter	Method	Specification	Result
Molecular weight		ca. 77 000	
Appearance		off-white to slightly pink lyophil.	corresponds
Protein content (%)	calc.on dried basis	min. 98.0	102
pH	3 % solution	6.5 – 8.0	7.15
Total bound Iron (µg/g)	ICP	max. 50	6.3
Endotoxin (EU/mg)	Limulus Amebocyte Lysate	max. 1.0	0.130
Moisture (%)	Loss on drying	max. 5.0	0.67
Minimum shelf life			05.04.2022
Storage (°C)			-15 to -25

This product was heat treated at 62°C for 10 hours for virus inactivation prior to lyophilization. Each unit of plasma used in the manufacture of this product has been tested and found to be non-reactive for Hepatitis B Surface Antigen (HbsAg), antibody to HIV I, II and anti-HCV, and Syphilis by FDA approved tests.

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle. We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.06.2018

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.
 This report has been computer generated and does not contain a signature.