

CERTIFICATE OF ANALYSIS

Trypsin NB Premium Grade, MS approved from porcine pancreas	Cat. No. : 37284 Lot. No. : 160578
---	---

Parameter	Method	Specification	Result
Sequence specificity	Digest of Insulin B chain (ox.) HPLC analysis	Must comply	Complies
Application (MS) approval	In-gel digestion of BSA and Mass spectrometric analysis	Must comply	Complies

For customer information:

Origin	Porcine pancreas
Molecular weight	24 kDa (trypsin)
Appearance	White to off-white powder
Purity (HPLC)	> 90%
Tryptic activity	> 6.000 U/g lyophilisate
Minimum shelf life	11.06.2018
Storage conditions	-15 to -25 °C

Sequence specificity assay:

Detection of two distinct peaks in Reversed Phase HPLC after digestion of oxidized B-chain of insulin at 37 °C for 18 h.

Application (MS) approval:

In-gel digest and identification of BSA.

Unit definition:

Tryptic activity:

1 Unit catalyzes the hydrolysis of 1 mmol Na-Benzoyl-L-arginine -4-nitroanilide hydrochloride (BAPNA) per minute at 30 °C, pH 8.0.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 12.10.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.