

## CERTIFICATE OF ANALYSIS

<b>Trypsin, Peptide (TP) Standard</b> from porcine pancreas	<b>Cat. No. :</b> 37285 <b>Lot. No. :</b> 150258
--	---

Parameter	Method	Specification	Result
Sequence specificity	Digest of Insulin B chain (ox.) HPLC analysis	Must comply	Complies

### For customer information:

Origin	Porcine pancreas
Molecular weight	24 kDa (trypsin)
Appearance	White to off-white powder
Purity (HPLC)	> 90%
Tryptic activity	> 6.000 U/g lyophilisate
Minimum shelf life	09/ 2016
Storage conditions	-15 to -25 °C

### Specificity assay:

Detection of two distinct peaks in Reversed Phase HPLC after digestion of oxidized B-chain of insulin at 37°C for 18h.

### Unit definition:

#### Tryptic activity:

1 Unit catalyzes the hydrolysis of 1 µmol N $\alpha$ -Benzoyl-L-arginine -4-nitroanilide hydrochloride (BAPNA) per minute at 30 °C, pH 8.0.

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 10.03.2015

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.