

# CERTIFICATE OF ANALYSIS

<b>Trypsin MS approved</b> from porcine pancreas	<b>Cat. No. :</b>	<b>37286</b>
	<b>Lot. No. :</b>	<b>170741</b>

Parameter	Method	Specification	Result
<b>Sequence specificity</b>	Digest of Insulin B chain (ox.) HPLC analysis	Must comply	Complies
<b>Application (MS) approval</b>	In-gel digestion of BSA and Mass spectrometric analysis	Must comply	Complies

### For customer information:

<b>Origin</b>	Porcine pancreas
<b>Molecular weight</b>	24 kDa (trypsin)
<b>Appearance</b>	White to off-white powder
<b>Purity (HPLC)</b>	> 90%
<b>Tryptic activity</b>	> 6.000 U/g lyophilisate
<b>Minimum shelf life</b>	26.01.2019
<b>Storage conditions</b>	-15 to -25 °C

### Sequence specificity assay:

Detection of two distinct peaks in Reversed Phase HPLC after digestion of oxidized B-chain of insulin at 37 °C for 18 h.

### Application (MS) approval:

In-gel digest and identification of BSA.

### Unit definition:

#### Tryptic activity:

1 Unit catalyzes the hydrolysis of 1 mmol Na-Benzoyl-L-arginine -4-nitroanilide hydrochloride (BAPNA) per minute at 30 °C, pH 8.0.

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 08.08.2017

Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

This report has been computer generated and does not contain a signature.