

**CERTIFICATE OF ANALYSIS**

<b>Trypsin 1:250 from bovine pancreas</b>	<b>Cat.No. :</b>	<b>37289</b>
<b>lyophil.</b>	<b>Contr.No. :</b>	<b>171027</b>

Parameter	Method	Specification	Result
<b>Appearance</b>		white lyophilisate	corresponds
<b>Activities:</b>			
<b>Trypsin (U/mg)</b>	Bergmeyer		3.3
<b>Trypsin (NF/USP U/mg)</b>	USP	≥ 250	260
<b>Chymotrypsin (NF/USP U/mg)</b>	USP	≤ 75	40
<b>Microbiological analysis</b>	Total viable count	£ 10 000 cfu/g	corresponds
	Pseudomonas aeruginosa	negative	corresponds
	Salmonella species	negative	corresponds
	Staphylococcus aureus	negative	corresponds
	E.coli	negative	corresponds
<b>Application</b>		suitable for cell culture	
<b>Minimum shelf life</b>			26.09.2019
<b>Storage (°C)</b>			-15 bis/ to -25

**Unit definition**

**Trypsin (Bergmeyer)**

1 U catalyzes the hydrolysis of 1 µmol N- $\alpha$ -benzoyl-L-arginine ethyl ester (BAEE) per minute at 25°C, pH 8.0.

**We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.**

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 16.11.2017

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.