

CERTIFICATE OF ANALYSIS

Trypsin 1:250 from bovine pancreas lyophil.	Cat.No. : 37289
	Contr.No. : 180851

Parameter	Method	Specification	Result
Appearance		white lyophilisate	corresponds
Activities:			
Trypsin (U/mg)	Bergmeyer		6.3
Trypsin (NF/USP U/mg)	USP	≥ 250	395
Chymotrypsin (NF/USP U/mg)	USP	≤ 75	66
Microbiological analysis	Total viable count	≤ 10 000 cfu/g	corresponds
	Pseudomonas aeruginosa	negative	corresponds
	Salmonella species	negative	corresponds
	Staphylococcus aureus	negative	corresponds
	E.coli	negative	corresponds
Application		suitable for cell culture	
Minimum shelf life			24.06.2020
Storage (°C)			-15 bis/ to -25

Unit definition

Trypsin (Bergmeyer)

1 U catalyzes the hydrolysis of 1 µmol N-α-benzoyl-L-arginine ethyl ester (BAEE) per minute at 25°C, pH 8.0.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 09.08.2018

Christian Monsler

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.