

CERTIFICATE OF ANALYSIS

Trypsin inhibitor from soybean
lyophil.

Cat.No. : 37329
Contr.No.: 191077

Parameter	Method	Specification	Result
Molecular weight		ca. 22 000	
Appearance		pale yellow lyophilisate	corresponds
Solubility	1 mg/ml in 0.067 M potassium-phosphate buffer, pH 7.6	clear and colorless	corresponds
Activity (U/mg)		min. 13 000	14 300
Minimum shelf life			05.09.2022
Storage (°C)			+2 to +8

Unit definition

1 IU (inhibitor units) inhibits 1 U trypsin as defined by cleavage of 1 µmole BAEE (N-benzoyl-L-arginine ethyl ester) per minute.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 11.11.2019

Christian Monsler

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.