

CERTIFICATE OF ANALYSIS

Urease from Jack bean lyophil.	Cat.No.: 37799
	Contr.No.: 190572

Parameter	Method	Specification	Result
Molecular weight		ca. 545 000	
Appearance		white lyophilisate	corresponds
Activity (U/mg)	Urease by Nessler's	min. 220	302
Activity in other units (U/mg)	Urease by Bergmeyer	min. 90	126
Minimum shelf life			04.06.2021
Storage (°C)			+2 to +8

Unit definitions

Urease (by Nessler's)

1 Unit is that amount of enzyme causing the liberation of 1 µmol ammonia per minute at 25°C, pH 7.0.

Urease (by Bergmeyer)

1 Unit is that amount of enzyme causing the hydrolysis of 1 µmol of urea per minute at 25°C, pH 8.0.

We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Quality Control

Printing date: 06.06.2019

Christian Monsler

Marina Konter

This report has been computer generated and does not contain a signature.