

# CERTIFICATE OF ANALYSIS

<b>Urease from Jack bean</b> lyophil.	<b>Cat.No. :</b> 37799
	<b>Contr.No.:</b> 210606

Parameter	Method	Specification	Result
<b>Molecular weight</b>		ca. 545 000	
<b>Appearance</b>		white lyophilisate	corresponds
<b>Activity (U/mg)</b>	Urease by Nessler's	min. 220	365
<b>Activity in other units (U/mg)</b>	Urease by Bergmeyer	min. 90	152
<b>Minimum shelf life</b>			17.05.2023
<b>Storage (°C)</b>			+2 to +8

## Unit definitions

### Urease (by Nessler's)

1 Unit is that amount of enzyme causing the liberation of 1  $\mu\text{mol}$  ammonia per minute at 25 °C, pH 7.0.

### Urease (by Bergmeyer)

1 Unit is that amount of enzyme causing the hydrolysis of 1  $\mu\text{mol}$  of urea per minute at 25 °C, pH 8.0.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 19.05.2021

Christian Monsler

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.