

**CERTIFICATE OF ANALYSIS**

---

<b>Phosphatase Inhibitor Mix I, powder</b>	<b>Cat.No. : 39050</b>
	<b>Lot No.: P090143</b>

---

<b>Components</b>	<b>Target-Phosphatase</b>
<b>Imidazol</b>	Alkaline Phosphatase
<b>Sodium fluoride</b>	Acid Phosphatase Protein-Phosphatase-2A Protein-Phosphatase-2B Protein-Phosphatase-2C
<b>Sodium molybdate</b>	Acid Phosphatase Phosphoprotein Phosphatase
<b>Sodium ortho vanadate</b>	Protein-Tyrosine-Phosphatase Alkaline Phosphatase
<b>Sodium tartrate</b>	Acid Phosphatase
<b>Storage (°C)</b>	+15 to +30
<b>Minimum shelf life</b>	05/ 2014
<b>Reconstitution</b>	Dissolve the content of the vial (156 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for treatment of 100 ml cell extract). After reconstitution, store at +2°C to +8°C (stable for 2 years).

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

---

**SERVA Electrophoresis GmbH**  
**Quality Control****Printing date:** 21.03.2013

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.