

CERTIFICATE OF ANALYSIS

Phosphatase Inhibitor Mix I, powder

Cat.No. : 39050

Lot No.: P130337

Components	Target-Phosphatase
Imidazol	Alkaline Phosphatase
Sodium fluoride	Acid Phosphatase Protein-Phosphatase-2A Protein-Phosphatase-2B Protein-Phosphatase-2C
Sodium molybdate	Acid Phosphatase Phosphoprotein Phosphatase
Sodium ortho vanadate	Protein-Tyrosine-Phosphatase Alkaline Phosphatase
Sodium tartrate	Acid Phosphatase
Storage (°C)	+15 to +30
Minimum shelf life	13.08.2017
Reconstitution	Dissolve the content of the vial (156 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for treatment of 100 ml cell extract). After reconstitution, store at +2°C to +8°C (stable for 2 years).

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 17.05.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.