

# CERTIFICATE OF ANALYSIS

**Phosphatase Inhibitor Mix I, powder**

**Cat.No. : 39050**

**Lot No.: P220068**

<b>Components</b>	<b>Target-Phosphatase</b>
<b>Imidazol</b>	Alkaline Phosphatase
<b>Sodium fluoride</b>	Acid Phosphatase Protein-Phosphatase-2A Protein-Phosphatase-2B Protein-Phosphatase-2C
<b>Sodium molybdate</b>	Acid Phosphatase Phosphoprotein Phosphatase
<b>Sodium ortho vanadate</b>	Protein-Tyrosine-Phosphatase Alkaline Phosphatase
<b>Sodium tartrate</b>	Acid Phosphatase
<b>Storage (°C)</b>	+15 to +30
<b>Minimum shelf life</b>	28.02.2026
<b>Reconstitution</b>	Dissolve the content of the vial (156 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for treatment of 100 ml cell extract). After reconstitution, store at +2°C to +8°C (stable for 2 years).

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**

**Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 14.05.2024**

Dr. Judith Koch

Daniela Lux-Helmstetter

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**

This report has been computer generated and does not contain a signature.