

CERTIFICATE OF ANALYSIS

Phosphatase Inhibitor Mix I, powder

Cat.No. : 39050

Lot No.: P250337

Components	Target-Phosphatase
Imidazol	Alkaline Phosphatase
Sodium fluoride	Acid Phosphatase Protein-Phosphatase-2A Protein-Phosphatase-2B Protein-Phosphatase-2C
Sodium molybdate	Acid Phosphatase Phosphoprotein Phosphatase
Sodium ortho vanadate	Protein-Tyrosine-Phosphatase Alkaline Phosphatase
Sodium tartrate	Acid Phosphatase

Storage (°C) +15 to +30

Minimum shelf life 15.10.2027

Reconstitution Dissolve the content of the vial (156 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for treatment of 100 ml cell extract). After reconstitution, store at +2°C to +8°C (stable for 2 years).

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 21.10.2025

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.