

# CERTIFICATE OF ANALYSIS

**Phosphatase Inhibitor Mix II**  
solution

**Cat.No. : 39055**  
**Lot No.: P120483**

<b>Components</b>	<b>Target-Phosphatase</b>
<b>EDTA</b>	Alkaline Phosphatase Protein-Phosphatase-2B Protein-Phosphatase-2C
<b>di-Sodium-<math>\beta</math>-Glycerophosphate</b>	Serine/Threonine- Phosphatase
<b>Imidazol</b>	Alkaline Phosphatase
<b>Sodium fluoride</b>	Acid Phosphatase Protein-Phosphatase-2A Protein-Phosphatase-2B Protein-Phosphatase-2C
<b>Sodium molybdate</b>	Acid Phosphatase Phosphoprotein Phosphatase
<b>Sodium ortho vanadate</b>	Protein-Tyrosine-Phosphatase Alkaline Phosphatase
<b>Sodium tartrate</b>	Acid Phosphatase
<b>Storage (°C)</b>	+2°C to +8°C
<b>Minimum shelf life</b>	11/ 2014

Aqueous stock solution (100x): 1 ml is suitable  
for treatment of 100 ml tissue extract.

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date: 03.12.2012**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

This report has been computer generated and does not contain a signature.