

CERTIFICATE OF ANALYSIS

Phosphatase Inhibitor Mix II
solution

Cat.No. : 39055
Lot No.: P130318

Components	Target-Phosphatase
EDTA	Alkaline Phosphatase Protein-Phosphatase-2B Protein-Phosphatase-2C
di-Sodium-β-Glycerophosphate	Serine/Threonine- Phosphatase
Imidazol	Alkaline Phosphatase
Sodium fluoride	Acid Phosphatase Protein-Phosphatase-2A Protein-Phosphatase-2B Protein-Phosphatase-2C
Sodium molybdate	Acid Phosphatase Phosphoprotein Phosphatase
Sodium ortho vanadate	Protein-Tyrosine-Phosphatase Alkaline Phosphatase
Sodium tartrate	Acid Phosphatase
Storage (°C)	+2°C to +8°C
Minimum shelf life	01.08.2017

Aqueous stock solution (100x): 1 ml is suitable
for treatment of 100 ml tissue extract.

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 17.05.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.