

CERTIFICATE OF ANALYSIS

Phosphatase Inhibitor Mix II
solution

Cat.No. : 39055
Lot No.: P230402

Components	Target-Phosphatase
EDTA	Alkaline Phosphatase Protein-Phosphatase-2B Protein-Phosphatase-2C
di-Sodium- β -Glycerophosphate	Serine/Threonine- Phosphatase
Imidazol	Alkaline Phosphatase
Sodium fluoride	Acid Phosphatase Protein-Phosphatase-2A Protein-Phosphatase-2B Protein-Phosphatase-2C
Sodium molybdate	Acid Phosphatase Phosphoprotein Phosphatase
Sodium ortho vanadate	Protein-Tyrosine-Phosphatase Alkaline Phosphatase
Sodium tartrate	Acid Phosphatase
Storage (°C)	+2°C to +8°C
Minimum shelf life	09.11.2025

Aqueous stock solution (100x): 1 ml is suitable
for treatment of 100 ml tissue extract.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 14.11.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.