

**CERTIFICATE OF ANALYSIS**

<b>Protein Molecular Weight Standards</b> range 0.2 – 450 kDa	<b>Kat.Nr./Cat.No. :</b>	<b>39064</b> <b>P060155</b>
--	--------------------------	--------------------------------

<b>Protein Components</b>	<b>Molecular weight</b>
DNP-L-Alanine	255
Cytochrome C pig heart	12 400
Myoglobin equine	17 800
Chymotrypsinogen A	25 000
Albumin egg	45 000
Albumin bovine (BSA)	67 000
Aldolase rabbit	160 000
Catalase bovine	240 000
Ferritin horse	450 000
Storage (°C)	-15 to -25
Minimum shelf life	31.12.2007

**The kit MS II contains 25 mg of each of the above mentioned, water soluble substances with well defined molecular weights in native electrophoresis.**

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.**  
We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Qualitätskontrolle/ Quality Control**

**Druckdatum/ printing date: 01.04.2016**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

**Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.**  
This report has been computer generated and does not contain a signature.