

CERTIFICATE OF ANALYSIS

Protein Molecular Weight Standards range 12.4 – 450 kDa for native PAGE	Cat. No. :	39064
	Lot. No. :	P150109

Protein Components	Molecular weight
Cytochrome C horse heart	12 400
Myoglobin equine	17 800
Chymotrypsinogen A	25 000
Albumin egg	45 000
Albumin bovine (BSA)	67 000
Aldolase rabbit	160 000
Catalase bovine	240 000
Ferritin horse	450 000
Storage (°C)	-15 to -25
Minimum shelf life	22.11.2016

The kit MS II contains 25 mg of each of the above mentioned, water soluble substances with well defined molecular weights in native electrophoresis.

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 10.11.2015

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.