

**CERTIFICATE OF ANALYSIS**

<b>Protein Molecular Weight Standards</b> range 12.4 – 450 kDa for native PAGE	<b>Cat. No. :</b>	<b>39064</b>
	<b>Lot. No. :</b>	<b>P150438</b>

<b>Protein Components</b>	<b>Molecular weight</b>
<b>Cytochrome C horse heart</b>	12 400
<b>Myoglobin equine</b>	17 800
<b>Chymotrypsinogen A</b>	25 000
<b>Albumin egg</b>	45 000
<b>Albumin bovine (BSA)</b>	67 000
<b>Aldolase rabbit</b>	160 000
<b>Catalase bovine</b>	240 000
<b>Ferritin horse</b>	450 000
<b>Storage (°C)</b>	-15 to -25
<b>Minimum shelf life</b>	22.11.2016

**The kit MS II contains 25 mg of each of the above mentioned, water soluble substances with well defined molecular weights in native electrophoresis.**

**We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.**

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date: 19.11.2015**

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.