

CERTIFICATE OF ANALYSIS

Protein Molecular Weight Standards	Cat. No. :	39064
range 12.4 – 450 kDa for native PAGE	Lot. No. :	P170130

Protein Components	Molecular weight
Cytochrome C horse heart	12 400
Myoglobin equine	17 800
Chymotrypsinogen A	25 000
Albumin egg	45 000
Albumin bovine (BSA)	67 000
Aldolase rabbit	160 000
Catalase bovine	240 000
Ferritin horse	450 000
Storage (°C)	-15 to -25
Minimum shelf life	30.07.2018

The kit MS II contains 25 mg of each of the above mentioned, water soluble substances with well defined molecular weights in native electrophoresis.

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 13.12.2017

Dipl.-Ing. (FH) Bernhard Göckel

Marina Konter

This report has been computer generated and does not contain a signature.