

# CERTIFICATE OF ANALYSIS

<b>Protein Molecular Weight Standards</b> range 12.3 – 450 kDa for native PAGE	<b>Cat. No. :</b> 39064 <b>Lot. No. :</b> P230426
---	--

<b>Protein Components</b>	<b>Molecular weight</b>
<b>Cytochrome C</b> (porcine, heart)	12 300
<b>Myoglobin</b> (equine)	17 800
<b>Chymotrypsinogen A</b>	25 000
<b>Albumin</b> (egg)	45 000
<b>Albumin</b> (bovine, BSA)	67 000
<b>Aldolase</b> (rabbit)	160 000
<b>Catalase</b> (bovine)	240 000
<b>Ferritin</b> (horse)	450 000
<b>Storage (°C)</b>	-15 to -25
<b>Minimum shelf life</b>	29.10.2024

The kit contains 25 mg of each of the above mentioned, water soluble substances with well defined molecular weights in native electrophoresis.

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 12.12.2023

Dr. Judith Koch

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.