

# CERTIFICATE OF ANALYSIS

**Protease-Inhibitor-Mix G**

**Cat.No. : 39101**

**Lot No.: P210149**

<b>Cat. number</b>	<b>Components</b>	<b>Target Protease</b>
12745	AEBSF-HCl	Serine Proteases
13718	Aprotinin from bovine lung	Serine Proteases
21100	E-64	Cysteine Proteases
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
11280	EDTA-disodium	Metalloproteases
<b>Storage (°C)</b>		-15 to -25
<b>Minimum shelf life</b>		04.05.2024
<b>Reconstitution</b>	Dissolve the content of the vial (31.2mg) in 1ml distilled water to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for 4 weeks).	

We do not guarantee that the product can be used for a special application.  
This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH**  
**Quality Control**

**Printing date:** 04.05.2021

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.