

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix M	Cat.No. : 39102
	Lot No.: P160245

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
13718	Aprotinin from bovine lung	Serine Proteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
21100	E-64	Cysteine Proteases
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
52682	Pepstatin A	Aspartic Proteases

Storage (°C) -15 to -25

Minimum shelf life 14.06.2019

Reconstitution Dissolve the content of the vial (15.8mg) in 1ml DMSO to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).

**We do not guarantee that the product can be used for a special application.
This document does not release you from performing the standard control upon receipt of incoming goods.**

SERVA Electrophoresis GmbH
Quality Control

Printing date: 14.06.2016

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.