

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix M

Cat.No. : 39102

Lot No.: P200316

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
13718	Aprotinin from bovine lung	Serine Proteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
21100	E-64	Cysteine Proteases
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
52682	Pepstatin A	Aspartic Proteases
Storage (°C)		-15 to -25
Minimum shelf life		09.09.2023

Reconstitution Dissolve the content of the vial (15.8mg) in 1ml DMSO to obtain a 100x stock solution (suitable for the treatment of 100ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.

Dieses Dokument entbindet nicht von einer branchenüblichen Eingangskontrolle.

We do not guarantee that the product can be used for a special application.

This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH
Qualitätskontrolle/ Quality Control

Druckdatum/ printing date: 09.09.2020

Christian Monsler

Patricia Milford

Dieses Dokument ist ein Computerausdruck und deshalb ohne Unterschrift gültig.

This report has been computer generated and does not contain a signature.