

CERTIFICATE OF ANALYSIS

Protease-Inhibitor-Mix P	Cat.No. : 39103
	Lot No.: P130055

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
21100	E-64	Cysteine Proteases
52682	Pepstatin A	Aspartic Proteases
32010	1,10-Phenanthroline	Metalloproteases
14980	Bestatin	Aminopeptidase B and Leucine Aminopeptidase
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
Storage (°C)		-15 to -25
Minimum shelf life		02/ 2016

Reconstitution Dissolve the content of the vial (64.2 mg) in 1ml DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).

<p>We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.</p>
--

SERVA Electrophoresis GmbH
Quality Control

Printing date: 20.02.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann

This report has been computer generated and does not contain a signature.